

PRODUCT LISTING

Classification	REF No.	Name	Specs	
Nucleic Acid Electrophoresis	CP17201M	TAE	100 pouches	
	CP17201-50 L	TAE	1 bottle	
	CP17202M	TBE	50 pouches	
	CP23201M	Tris-EDTA (pH 8.0)	100 pouches	
Protein Electrophoresis	CP17204M	Tris-MES-SDS	50 pouches	
	CP17205M	Tris-MOPS-SDS	50 pouches	
	CP17206M	Tris-Glycine-SDS	50 pouches	
	CP17206-50 L	Tris-Glycine-SDS	1 bottle	
	CP21204M	Tris-HEPES-SDS	100×0.5 L	
	CP23301S	30% Acr-Bis (29:1)	500 ml	
	Immunoblotting	CP17207M	PBS-T	100 pouches
CP17207-50 L		PBS-T	1 bottle	
CP17208M		PBS	100×1 L	
CP17208-50 L		PBS	1 bottle	
CP17209M		TBS-T	100 pouches	
CP17209-50 L		TBS-T	1 bottle	
CP17210M		TBS	100 pouches	
CP17210-50 L		TBS	1 bottle	
CP20203L		Rapid Blocking Buffer (TBS-T)	20 pouches	
CP20204M		Transfer Buffer(Semi Dry)	100 pouches	
CP21202M		Tris-Glycine	50 pouches	
CP21211M		PBS	50×2 L	
Dust-free Granules		CP21209S	SDS dust-free Granules	250 g
		CP21209M	SDS dust-free Granules	1000 g



Welcome to negotiate custom specifications or formulations.



Scan and explore for more!

BestEnzymes Biotech Co.,Ltd.

Add: Building 1, Floors 12-14, 17 Huaguo Mountain Avenue, Gaoxin Technology Industrial Development Zone, Haizhou District, Lianyungang City, Jiangsu Province.
Tel: +86-18861343080 · Tina@best-enzymes.com · <http://en.www.best-enzymes.com>

ver. 202401



Buffer Solution Instant Granules

Strong Homogeneity/Rapid Dissolution/
Reliable Quality/Ready to Use



Best Enzymes For Better Life

<http://en.www.best-enzymes.com>


Advantages

• Stronger homogeneity


Through the process of granulation, the powder is thoroughly mixed and formed into small instant granules that have consistent composition. In contrast, conventional powders may have variations in component composition due to differences in particle size and density. Therefore, instant granules have stronger homogeneity.

• Faster dissolution rate

When powdered substances are dissolved, they tend to clump together, forming a high-concentration liquid layer that encapsulates the powder granules, impeding their dissolution. However, after being processed into instant granules, each granule contains numerous pores or channels, leading to an increased surface area and significantly enhancing the dissolution rate. This improved surface area-to-volume ratio allows for faster and more efficient dissolution of the instant dissolve granules.



Selection in different scenarios of nucleic acid electrophoresis.




Electrophoresis conditions: TAE/TBE—6 V/cm, 40 min; Rapid Running Buffer—15 V/cm, 15 min

- It is considered that TAE buffer provides better separation for large DNA fragments, while TBE buffer is more effective in separating small DNA fragments.
- Using Rapid Running Buffer can save over 50% of electrophoresis time.

Buffer Solution Instant Granules

Instant granules for all processes of Western Blot



Protein Electrophoresis Buffer Selection

Molecular Weight (kDa)	8%			10%			12%			4~12%		
	MES	MOPS	HEPES	MES	MOPS	HEPES	MES	MOPS	HEPES	MES	MOPS	HEPES
200			200			200			200			200
116		200		200			116			200		
97.4			116	116		116	97.4		116	97.4		116
66.3			97.4	66.3	97.4	97.4	66.3		66.3	116		97.4
55.4		116		55.4			55.4		55.4	97.4		
			66.3			66.3			66.3			66.3
		66.3		36.5	55.4	55.4	36.5		55.4	66.3		66.3
			55.4			36.5			36.5	36.5		55.4
		55.4		21.5		31			31	21.5		55.4
			36.5			31			31	14.4		36.5
				14.4		6			14.4	6		31
						21.5			21.5	14.4		36.5
						6			6	31		31
										6		21.5
										3.5		14.4
										2.5		14.4
												14.4

① Bis-Tris Gel System

- The MES buffer is suitable for low to medium molecular weight proteins (6~260 kDa).
- The MOPS buffer is suitable for medium to high molecular weight proteins (14~260 kDa).
- The HEPES buffer exhibits more uniform migration and has a wide range of molecular weights.
- The HEPES buffer without SDS is also provided for native gels.

② Tris-Glycine Gel System

- The Tris-Glycine-SDS (TGS) buffer is suitable for proteins with a wide range of molecular weights (6~400 kDa).
- The TG buffer, which is a version of the TGS buffer without SDS, is suitable for native gels.