

# mRNA *In vitro* Synthesis Solutions

Template Generation • Template Linearization  
*In vitro* Transcription • RNA Analysis

---



# mRNA *In vitro* Synthesis Solutions


Products manufactured in compliance with GMP regulations can meet diverse customer requirements.

## GMP Raw Material Quality Standards

Quality Control	Quality Standard	Detection Method
Appearance	Colorless clear liquid	Visual Check
Protein Purity	≥ 95%	SDS-PAGE
Non-specific Endonuclease Activity	Supercoiled plasmid conformational change < 20%	20 U enzyme was incubated with supercoiled plasmid at 37°C for 4 hours
DNase Activity	No dsDNA degradation observed	20 U enzyme was incubated with dsDNA fragments at 37°C for 16 hours
RNase Activity	RNA degradation ≤ 10%	20 U enzyme was incubated with RNA at 37°C for 1 hour CFU/mL
Host Cell DNA	< 10 copies/20 U	ChP(2025) Volume IV, Detection of Exogenous DNA Residues, Method 3 (General Chapter 3407)
Host Cell Protein	< 50 ppm	ChP(2025) Volume IV, Measurement of E. coli Host Cell Proteins (General Chapter 3412)
Bacterial Endotoxin	< 10 EU/mg	ChP(2025) Volume IV, Bacterial Endotoxins Test, Method 1 (General Chapter 1143)
Mycoplasma	Negative	Mycoplasma test kit
Heavy Metals	< 10 ppm	ChP(2025) Volume IV, Limit Test for Heavy Metals, Method I (General Chapter 0821)
Microbial Limit	The total aerobic microbial count is below 5 cfu/ml, and the total combined yeasts/molds count is below 5 cfu/ml.	ChP(2025) Volume IV, Microbial Limit Tests for Non-Sterile Products: Microbial Enumeration Method (General Chapter 1105)

Note: The above information is given as examples of Bsal, GMP Grade .

## mRNA Vaccine Manufacturing Process



BestEnzymes Biotech offer a wide range of GMP and non-GMP grade raw materials, assisting in the research and production of mRNA vaccines.

## mRNA Synthesis Workflow Example & Available Products

1	2	3	4	5
Template Generation	Template Linearization	<i>In vitro</i> Transcription	RNA Capping*and Tailing	RNA Analysis
Restriction Endonucleases	Bsal, GMP Grade	T7 High Yield RNA Synthesis Kit	Cap1 analogue AG	RNase A
2× S705 HiFi Master Mix	XbaI, GMP Grade	T7 RNA Polymerase	Poly(A) RNA Polymerase	RNase T1
Ultra / T4 DNA Ligase	XhoI, GMP Grade	High T7 RNA Polymerase		Thermostable RNase T1
DNA Assembly Mix Ultra/ Mono	NheI, GMP Grade	RNase Inhibitor, Murine		RNase H
Golden Gate Assembly Kit (BsmBI)	BspQI, GMP Grade	DNase I, RNase-free		Thermostable RNase H
Golden Gate Assembly Kit (BpiI)	Esp3I(BsmBI), ADCF	DNase I-ST, RNase-free		MazF
2× Taq PCR Aurora Mix	EcoRI, ADCF	Pyrophosphatase, Inorganic (Yeast)		RNase GG
phi29 II DNA polymerase	Bsal ELISA Kit	Pseudo-UTP		
TeIn Protelomerase	XhoI ELISA Kit	N <sup>1</sup> -Methyl-Pseudo-UTP		
dNTPS	XbaI ELISA Kit	NTPs		
	BsmBI ELISA Kit			

\*Accurately controlling poly(A) tail length using enzymatic methods is challenging, which is why the poly(A) sequence is typically incorporated directly into the plasmid template.




# In vitro Transcription

The T7 High Yield RNA Synthesis Kit is capable of synthesizing long or short RNA transcripts in large quantities from DNA templates containing a T7 promoter. The kit's T7 Enzyme Mix is optimized with RNase inhibitor and inorganic pyrophosphatase, ensuring highly efficient and robust RNA synthesis. The synthesized RNA is applicable for a wide range of downstream applications, including *in vitro* translation, RNA structure and function studies, RNase protection assays, probe hybridization, RNA interference, and more.


## Product Feature

- High yield, 1 µg of linearized double-stranded DNA template can yield at least 150 µg of RNA.
- Fast rate, 80 minutes will be sufficient to obtain the required amount of the product.
- Wide adaptability to template lengths (100 bp to 9 kb), suitable for preparing mRNA, CRISPR-required sgRNA, and other applications.
- Optimized reaction system, compatible with modified NTPs and cap analogs.


## Performance Demonstration




The relationship between amount of template and product.



The relationship between reaction time of different templates length and product output.



The relationship between different template lengths and the products output.



As tested by customers, when using cap analogues for co-transcription, the integrity of the capped mRNA exceeds 85%.

# Template Generation

When amplifying the target fragment, it is essential to select a high-fidelity DNA polymerase to ensure the accuracy of the template. BestEnzymes Biotech has launched the 2×S705 HiFi Master Mix, a high-fidelity PCR premix containing S705 High-Fidelity DNA Polymerase. This polymerase offers 70-fold higher fidelity than Taq polymerase and 5-fold faster amplification than Pfu polymerase, meeting the requirements for accurate template amplification.

BestEnzymes Biotech offers a comprehensive catalog including over a hundred restriction endonucleases, thermostable reverse transcriptase, T4 DNA Ligase (Fast), and Gibson assembly cloning kits, effectively addressing the needs for plasmid template construction.

## Template Linearization

When utilizing plasmid DNA as an *in vitro* transcription template, it's crucial to select appropriate restriction endonucleases to ensure complete digestion and linearization of the circular plasmid. BestEnzymes Biotech offers a diverse range of GMP-grade enzymes. We guarantee full traceability throughout our entire production process, including all raw materials and excipients. Our enzymes are completely free of antibiotics and any animal-derived components. Furthermore, we rigorously control process-related impurities such as host proteins, exogenous DNA, non-specific endonucleases, DNase, and RNase, as well as microbial limits and bacterial endotoxins. This stringent quality control meets the demanding raw material and excipient requirements for critical applications like vaccine and drug production and cell therapy.

Meanwhile, BestEnzymes Biotech provides GMP-grade restriction endonucleases customization to meet users' personalized needs for plasmid templates.

To address the issue of restriction endonuclease residue after plasmid template linearization, BestEnzymes Biotech pioneered the industry's first series of detection kits for restriction endonuclease residues. Currently, ELISA kits for individual enzymes including BsaI, XhoI, XbaI, and BsmBI are available, with new products to be gradually launched according to customer demands.

## RNA Analysis

When performing *in vitro* transcription with both capping and tailing to generate complete mRNA, it's crucial to analyze the quality of the RNA product. To address this need, BestEnzymes Biotech offers a diverse array of RNase. These enzymes provide versatile options for RNA secondary structure analysis, capping rate evaluation, and comprehensive integrity assessment.


Product Name	Reaction Substrate	Recognition Sequence	Application
<b>RNase A</b>	Primarily single-stranded RNA	C/U	Preparation of plasmids and genomic DNA Remove RNA from recombinant protein preparations RNA In combination with RNase T1 Identify single-base mutations in DNA or RNA
<b>RNase T1/ Thermostable RNase T1</b>	Single-stranded RNA	G	RNA sequencing Analysis of RNA secondary structure Combined with RNase A to remove RNA mRNA quality analysis
<b>RNase H/ Thermostable RNase H</b>	The RNA chain in the RNA-DNA hybrid molecule	Positional preference	Analysis of mRNA cap addition efficiency Remove mRNA before synthesizing the second strand of cDNA Removal of the poly(A) sequence from mRNA after oligo(dT) hybridization
<b>MazF</b>	Single-stranded RNA	Unmodified ACA	mRNA analysis

## Ordering Information

Classification	REF No.	Name	Specs
<b>Template Construction</b>	EG15124S	M-MLV III Reverse Transcriptase	10000 U
	EG24110S	2× S705 HiFi Master Mix	1 ml
	EG24110M	2× S705 HiFi Master Mix	5×1 ml
	EG24110L	2× S705 HiFi Master Mix	20×1 ml
	EG15205S	T4 DNA Ligase (Fast)	1000 U
	EG25208S	Ultra T4 DNA Ligase	1000 U
	EG24204S	DNA Assembly Mix Ultra	50 rxns
	EG25202S/M	DNA Assembly Mix Mono	50 rxns/250 rxns
	EG25202M	DNA Assembly Mix Mono	250 rxns
	EG25207S	Golden Gate Assembly Kit (BsmBI)	50 rxns
	EG25209S	Golden Gate Assembly Kit (Bpil)	50 rxns
	EG25104M/L	2× Taq PCR Aurora Mix	5×1 ml/100×1 ml
	EG25109S/M	phi29 II DNA Polymerase	500 U/2500 U
	EG25203S/M	TeIN Protelomerase	250 U/1000 U
	EG20907S	dNTPs (10 mM each)	1 ml
	EG20901S	dATP (100 mM)	1 ml
	EG20902S	dTTP (100 mM)	1 ml
	EG20903S	dCTP (100 mM)	1 ml
	EG20904S	dGTP (100 mM)	1 ml
	EG20905S	dUTP (100 mM)	1 ml

## Ordering Information

Classification	REF No.	Name	Specs
Template Generation	GMP501S	Bsal, GMP Grade	20000 U
	GMP502S	Xbal, GMP Grade	20000 U
	GMP503S	BspQI, GMP Grade	10000 U
	GMP504S	XhoI, GMP Grade	20000 U
	GMP505S	NheI, GMP Grade	20000 U
	EG21504S	Esp3I(BsmBI), ADCF	2000 U
	EG21504H	Esp3I(BsmBI), ADCF	2000 U
	EG22504S	EcoRI, ADCF	5000 U
	EG24601S	Bsal ELISA Kit	96 T
	EG24602S	XhoI ELISA Kit	96 T
	EG24603S	XbaI ELISA Kit	96 T
	EG24604S	BsmBI ELISA Kit	96 T
	In vitro Transcription	EG24104S	T7 High Yield RNA Synthesis Kit
GMP101S		T7 RNA Polymerase, GMP Grade	20000 U
GMP101M		T7 RNA Polymerase, GMP Grade	200000 U
EG20124S		T7 RNA Polymerase	5000 U
EG20124M		T7 RNA Polymerase	25000 U
EG20125S		High T7 RNA Polymerase	5000 U
EG20125M		High T7 RNA Polymerase	25000 U
EG21916S		ATP (100 mM)	1 ml
EG21917S		CTP (100 mM)	1 ml
EG21918S		GTP (100 mM)	1 ml
EG21919S		UTP (100 mM)	1 ml
EG21920S		NTP Bundle (25 mM each)	1 ml
EG21921S		N1-Methyl-Pseudo-UTP (100 mM)	100 µl
EG21922S		Pseudo-UTP (100 mM)	200 µl
GMP102S		Murine RNase Inhibitor, GMP Grade	40000 U
EG20002S		RNase Inhibitor, Murine	10000 U
GMP105S		DNase I-ST, GMP Grade	2000 U
EG24202S		DNase I-ST, RNase-free	1000 U
EG24206S		DNase I-ST, RNase-free (Lyophilized)	1000 U
GMP103S		Pyrophosphatase, Inorganic (Yeast, GMP Grade)	100 U
EG23109S	Pyrophosphatase, Inorganic (Yeast)	10 U	
EG23109M	Pyrophosphatase, Inorganic (Yeast)	100 U	
RNA Capping	EG22909S	Cap1 analogue AG	100 µl
	EG25108S	Poly(A) RNA Polymerase	100 U
	EG25108M	Poly(A) RNA Polymerase	500 U
RNA Analysis	EG20003S	RNase A	100 mg
	EG20003M	RNase A	1 g
	EG24210S	RNase T1	100 KU
	EG24210M	RNase T1	5×100 KU
	EG24211S	Thermostable RNase T1	100 KU
	EG24211M	Thermostable RNase T1	5×100 KU
	EG20208S	RNase H	1000 U
	EG24201S	Thermostable RNase H	200 U
	EG24505S	MazF	1000 U
	EG25521S	RNase GG	2500 U
EG25521M	RNase GG	12500 U	



WeChat Official Account



Research Subscription Account



WeChat Channels

## BestEnzymes Biotech Co., Ltd.

Add: No.17 Huaguoshan Avenue, Lianyungang City, Jiangsu Province, China  
Tel: 0518-8558 6628 · support@best-enzymes.com · <http://en.best-enzymes.com>

