

## TEV Protease (with His-tag)

REF: EG24303S

### Storage Condition

-20 °C , It is recommended to aliquot and store upon first use to avoid repeated freeze-thaw cycles.

### Components

Component	Amount
TEV Protease (10 U/μl)	100 μl
10× TEV Buffer	1 ml

### Description

TEV Protease (with His-tag) is derived from the tobacco etch virus and has been genetically engineered for recombinant expression in *Escherichia coli*. It retains the natural enzymatic activity of TEV protease while exhibiting improved tolerance to variations in pH and temperature. TEV Protease is a highly specific cysteine protease, widely utilized for the removal of affinity tags from fusion proteins. It precisely recognizes the seven-amino-acid sequence EXXYXQ↓(G/S) and cleaves between the glutamine and glycine/serine residues. The most commonly recognized sequence is Glu-Asn-Leu-Tyr-Phe-Gln↓Gly. The inclusion of a His-tag at the N-terminus of TEV Protease allows for its efficient removal via Ni-affinity chromatography following digestion, thereby simplifying the purification of the target protein.

### Definition of Activity Unit

One unit of activity is defined as the amount of enzyme required to cleave >85% of 3 μg of substrate within 1 hour at 30°C in 1× TEV Buffer.

### Storage Buffer

50 mM Tris-HCl; 250 mM NaCl; 10 mM DTT; 1 mM EDTA; 50% Glycerol; pH 7.5 @ 25°C.

### Quality Control Assays

#### Protein Purity

The enzyme is ≥95% pure as determined by SDS-PAGE analysis using Coomassie Blue staining.

#### Non-specific protease activity

A 50 μl reaction containing 3 μg of standard protein and 10 U of TEV Proteinase incubated for 1 hour at 30 °C results in no detectable degradation of the protein mixture as determined by SDS-PAGE.

### Protocol

1. Prepare the following reaction system in an EP tube:

Reagent	Amount
Substrate	20 μg
10× TEV Buffer	5 μl
TEV Protease (10 U/μl)	1 μl
ddH <sub>2</sub> O	Up to 50 μl

2. Incubate at 37°C for 1 hour or Incubate at 4°C overnight (about 16 hours) .

3. Add 16 μl of Protein Loading Buffer (SDS, 4×) (REF: EG21323S) to the above EP tube. Boil the sample for 10 minutes, and then take 10 μl for SDS-PAGE analysis.

4. After evaluating the initial results, you may optimize the cleavage reaction for your specific protein by optimizing the amount of TEV Protease, incubation temperature, or reaction time.

### Notice

1. Compatible with the following protease inhibitors: aprotinin, benzamidine, leupeptin, pepstatin, PMSF.

2. Compatible with 400 mM imidazole, 2 M urea, 0.2 M NaCl, 10 mM MgSO<sub>4</sub>, 10 mM MnCl<sub>2</sub>, 10 mM CaCl<sub>2</sub>, and up to 100 mM EDTA.

3. Inhibition occurs in the presence of ≥40% glycerol, ≥5 mM Zn<sup>2+</sup>, ≥1 mM Cu<sup>2+</sup>, and ≥10 mM Co<sup>2+</sup>.

4. The optimal reaction conditions for TEV are pH 7.0 and 30°C . However, it remains active within the range of pH6~9 and 4~37 °C . The reaction conditions can be selected according to the specific requirements of the target protein.