

## Tli DNA Polymerase (exo-)

REF: EG26106S

### Storage Condition

Store at -20°C for 2 years.

### Components

Component	Amount
Tli DNA Polymerase (exo-) (2 U/ $\mu$ l)	100 $\mu$ l
10 $\times$ Tli Buffer	1.5 ml
MgSO <sub>4</sub> (100 mM)	1 ml

### Description

Tli DNA Polymerase (exo-) is derived from a hyperthermophilic bacterium isolated from a submarine thermal vent. The enzyme is genetically modified and recombinantly expressed in *Escherichia coli* (*E. coli*), with its 3'→5' exonuclease activity removed. This polymerase features extreme thermostability: its half-life is 23 hours at 95 °C and 8 hours at 100 °C. It also possesses strong strand displacement activity. Its properties and operating protocols are fully consistent with Deep Vent<sup>®</sup> (exo-) DNA Polymerase. This product is well suited for whole genome amplification technologies such as MALBAC.

### Definition of Activity Unit

One unit is defined as the amount of enzyme that can incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 75°C.

### Quality Control Assays

#### Protein Purity

The enzyme is  $\geq$ 95% pure as determined by SDS-PAGE analysis using Coomassie Blue staining.

#### Endonuclease Activity

A 20  $\mu$ l reaction containing 200 ng of supercoiled plasmid and 2 U of Tli DNA Polymerase (exo-) incubated for 4 hours at 37°C results in <20% conversion to the nicked or linearized form as determined by agarose gel electrophoresis.

#### DNase Activity

A 20  $\mu$ l reaction containing 15 ng of dsDNA fragments and 2 U of Tli DNA Polymerase (exo-) incubated for 16 hours at 37 °C results in no detectable degradation of the dsDNA fragments as determined by agarose gel electrophoresis.

#### RNase Activity

A 10  $\mu$ l reaction containing 500 ng of RNA and 2 U of Tli DNA Polymerase (exo-) incubated for 1 hour at 37°C results in >90% of the substrate RNA remains intact as determined by agarose.

#### Single-stranded DNase Activity

A 20  $\mu$ l reaction containing 0.2 pmol of single-stranded DNA probe and 2 U of Tli DNA Polymerase (exo-) incubated for 30 minutes at 37°C results in less than 10% probe degradation as measured by real-time PCR.

#### Residual Host DNA

The residual *E. coli* host cell DNA in this product is less than 1 copy/2 U.

### Notice

- To increase amplification yield, add MgSO<sub>4</sub> to the reaction system to a final concentration of 1~2 mM. Please note that excessive additional Mg<sup>2+</sup> may cause non-specific amplification. Determine whether to supplement MgSO<sub>4</sub> according to your experimental requirements.
- For your safety and health, please wear a lab coat, disposable gloves and a mask when performing the operation.